

## HG CTNG Assay Kit - Compatible Transport Media Update Technical Bulletin – 1<sup>th</sup> March 2022

### HG CTNG Assay Kit -Compatible Transport Media Update

HiberGene are pleased to confirm validation of the following swabs for use with the HiberGene HGCTNGR assay for *Chlamydia trachomatis* and *Neisseria gonorrhoea*.

- 3ml Copan UTM + Copan Flocked Swab (Copan, Cat. No. 306C)
- 3ml Puritan UniTranz + Flocked Swab PurFlock Ultra (Puritan Medical Products Cat. No. UT-306)

Testing performed by the HiberGene R&D department confirms that this transport medium type delivers assay performance equivalent to that of previously validated transport media when utilised with the assay. An update version of the IFU can be found on the distributor portal.



---

If you have any questions or would like further information on supply of this product, please liaise with your relevant HiberGene contact.