

HG CTNG-19 AMF Update

Technical Bulletin – 24th November 2020

HG CTNG AMF Update

ATTENTION – Please note a new version of the CTNG .amf file is available for the HG SWIFT Instrument.

The updated HG CTNG amf (HG CTNG_rc2.amf) contains improved sample-calling parameters which further decrease the risk of false Positive results.

To ensure that you are using the latest version, please complete the following task:

- Please load the updated HG CTNG-19 Assay Method File (HG CTNG-19_rc2.amf) onto the HG Swift instrument.

All Assay Method Files are also available for download at <http://www.hibergene.com/distributors/distributor-portal/>



If you have any questions or would like further information on HiberGene tests, please liaise with your relevant HiberGene contact:

Chris Sloane e: chris.Sloane@hibergene.com

Aisling Healy e: aisling.healy@hibergene.com

Cathal McQuaid e: cathal.mcquaid@hibergene.com