

## HG CTNG Assay Kit - Compatible Transport Media Update Technical Bulletin – 13<sup>th</sup> August 2021



### HG CTNG Assay Kit -Compatible Transport Media Update

HiberGene are pleased to confirm validation of the 3ml Copan UTM + Copan Flocked Swab (Copan, Cat. No. 306C) for use with the HiberGene HGCTNGR assay for *Chlamydia trachomatis* and *Neisseria gonorrhoea*.

Testing performed by the HiberGene R&D department confirms that this transport medium type delivers assay performance equivalent to that of previously validated transport media when utilised with the assay.



---

If you have any questions or would like further information on supply of this product, please liaise with your relevant HiberGene contact.