

HG C. difficile Technical Bulletin – 19th December 2016


HG C. difficile

Our newest test, HG C. difficile (Product code: HGCDIFFR), is now CE marked and available.

From a stool sample, the simple protocol detects a portion of the *tcdA* gene which is conserved in both A+/B+ and A-/B+ toxigenic strains of *C. difficile*. The test will give you a positive result in less than 60 minutes.

This kit contains: swabs, diluent buffer #1, diluent buffer #2 and reaction strips.

Attention: a software update is required for your HG SWIFT instrument to run this test. Please see our software update bulletin for instructions on how to install the update.



If you have any queries or questions please contact your regional product manager:
roisin.fry@hibergene.com - Asia Pacific/Canada/USA/Latin America
tristan.dent@hibergene.com - Europe/Africa/Middle East/India/Pakistan